

Custom Machined Molded Ptfе Teflon Parts Manufacturer For Laboratory Ito Fto Conductive Glass Cleaning Flower Basket

Item Number: PL-1013



Introduction

High-purity PTFE flower baskets for semiconductor & lab use. Chemical-resistant, custom designs available. Ideal for silicon wafers & glass substrates.

[Learn More](#)

Diameter*groove width*groove length*applicable capacity

60mm*1.5mm*16mm*250ml beaker/20 pcs

60mm*2.5mm*16mm*250ml beaker/20 pcs

60mm*1.5mm*21mm*250ml beaker/20 pcs

60mm*2.5mm*21mm*250ml beaker/20 pcs

80mm*1.5mm*26mm*500ml beaker/20 pcs

80mm*2.5mm*26mm*500ml beaker/20 pcs

80mm*1.5mm*31mm*500ml beaker/20 pcs

80mm*2.5mm*31mm*500ml beaker/20 pcs

90mm*2.5mm*41mm*1000ml beaker/22 pcs