

Custom Ptfе Measuring Cylinders For Advanced Scientific And Industrial Applications

Item Number: PL-1009



Introduction

Discover the superior performance of our Polytetrafluoroethylene (PTFE) measuring cylinders, engineered for precision and resilience in demanding environments.

[Learn More](#)

Capacity	Total Height	Outer Diameter	Weight
10ml	108mm	17mm	42g
20ml	108mm	24mm	64g
50ml	118mm	34mm	128g
100ml	138mm	40mm	158g
250ml	147mm	60mm	374g
500ml	160mm	79mm	460g
1000ml	180mm	100mm	699g
2000ml	/	/	/