

Custom Ptfе Parts Manufacturer For Teflon Containers And Components

Item Number: PT-1011



Introduction

High-precision PTFE containers for labs & industry. Chemical-resistant, customizable sizes. Ideal for semiconductor, medical & lab applications. Get a quote!

[Learn More](#)

Round Vessel (Outer Dia. * Inner Dia. * Inner Depth)	Square/Rectangular Dish (Inner L x W * Inner Depth)
70mm*60mm*15mm	60*60mm *15mm
50mm*40mm*15mm	80*80mm*15mm
60mm*50mm*15mm	90*90mm*15mm
80mm*70mm*15mm	100*100mm*15mm
90mm*80mm*15mm	110*110mm*15mm
100mm*90mm*15mm	120*120mm*15mm
110mm*100mm*15mm	
120mm*110mm*15mm	
130mm*120mm*15mm	
160mm*150mm*15mm	200*200*120mm (approx. 5L PTFE cylinder)
92*83*19mm	10L PTFE cylinder
102*93*20mm	15L PTFE cylinder