

Custom Ptfе Shallow Evaporating Dishes For Diverse Applications

Item Number: PL-1002



Introduction

High-purity PTFE shallow evaporating dish for labs. Chemical-resistant, non-stick, ideal for sample evaporation. Custom sizes available. Order now!

[Learn More](#)

Specification	Bottom Outer Diameter (mm)	Mouth Inner Diameter (mm)	Total Height (mm)	Bottom Thickness (mm)	Wall Thickness (mm)
Shallow Type 20mm	20	15	10	5	2.5
Shallow Type 30mm	30	25	10	5	2.5
Shallow Type 40mm	40	30	10	5	5
Shallow Type 50mm	50	40	10	5	5
Shallow Type 60mm	60	50	10	5	5
Shallow Type 70mm	70	60	10	5	5
Shallow Type 80mm	80	70	10	5	5
Shallow Type 90mm	90	80	10	5	5
Shallow Type 100mm	100	90	10	5	5