

Custom Ptfе Square Trays For Industrial And Laboratory Use

Item Number: PL-1005



Introduction

Premium PTFE square trays for labs & industry. Chemical-resistant, non-stick, high-purity PTFE trays for semiconductors & medical use. Custom sizes available.

[Learn More](#)

Model	Internal L*W*H (mm)	External L*W*H (mm)
40mm	40*40*20mm	30*30*15mm
60mm	60*60*20mm	50*50*15mm
80mm	80*80*20mm	70*70*15mm
100mm	100*100*20mm	90*90*15mm
120mm	120*120*20mm	110*110*15mm