

Custom Ptfе Teflon Parts Manufacturer Conductive Glass Substrate Cleaning Rack

Item Number: PL-1016



Introduction

High-purity PTFE flower baskets for lab cleaning, resistant to chemicals & extreme temps. Custom designs available for semiconductor & medical applications.

[Learn More](#)

10×10mm (slot e.g., 1.1mm)	20×20mm (slot e.g., 1.1mm)	25×25mm (slot e.g., 2.2mm)	30×30mm (slot e.g., 2.2mm)
15×15mm (slot e.g., 1.1mm)	25×25mm (slot e.g., 1.1mm)	30×30mm (slot e.g., 1.1mm)	