

Custom Ptfе Teflon Parts Manufacturer Ptfе Cleaning Rack

Item Number: PL-1015



Introduction

High-purity PTFE flower baskets for labs & semiconductor use. Chemical-resistant, -180°C to +250°C, custom sizes available. Contact KINTEK today!

[Learn More](#)

Standard Options: Diameter*Slot Width*Applicable Capacity

45*1.5mm*100ml

45*2.5mm*100ml

60*1.5mm*250ml

60*2.5mm*250ml

80*1.5mm*500ml

80*2.5mm*500ml

90*1.5mm*1000ml

90*2.5mm*1000ml