

Custom Ptfе Volumetric Flasks For Advanced Scientific And Industrial Use

Item Number: PL-1018



Introduction

Premium PTFE volumetric flasks for labs – chemical-resistant, non-stick, precise measurements. Ideal for semiconductor, medical & analytical applications. Shop now!

[Learn More](#)

Volume (mL)	Inner Mouth Diameter (mm)	Max Outer Diameter (mm)	Total Height (mm)	Weight (g)
25	17	41	97	41
50	23	52	117	73
100	23	65	117	111
200	29	88	158	232
250	29	90	169	277
500	34	110	213	409
1000	38	150	260	965