

Customizable Ptfе Crucibles For Laboratory And Industrial Applications

Item Number: PL-1006



Introduction

High-purity PTFE crucibles with lids for labs. Chemical-resistant, non-stick, and durable. Ideal for sample digestion, AA, ICP-MS. Custom sizes available.

[Learn More](#)

Capacity	Outer Top Diameter (mm)	Bottom Diameter (mm)	Total Height (mm)
20ml	38	30	55
25ml	39	30	55
30ml	48	30	56
50ml	51	33	67
100ml	61	43	70
250ml	79	50	107