

Customizable Ptfе Three Neck Flasks For Advanced Chemical Applications

Item Number: PL-1007



Introduction

Our high-performance Three-Neck Round-Bottom Flasks are expertly crafted from pure Polytetrafluoroethylene (PTFE), a material renowned for its exceptional properties, making it ideal for demanding laboratory and industrial environments.

[Learn More](#)

Capacity / Neck Configuration	Overall Height (approx.)
50ml	/
100ml	114mm
250ml/24*3	135mm
500ml/24*3	180mm
1000ml/24*3	/
2000ml/24*3	/