

Ptfe Deep Evaporating Dishes Customizable Laboratory And Industrial Solutions

Item Number: PL-1001



Introduction

High-quality PTFE deep evaporating dishes for labs. Chemical-resistant, non-stick, and customizable. Ideal for sample handling. Get yours now!

[Learn More](#)

Dish Type	Bottom Outer Diameter (mm)	Inner diameter of mouth (mm)	Height (mm)	Bottom Thickness (mm)	Wall Thickness (mm)
Deep Type 20mm	20	15	20	5	2.5
Deep Type 30mm	30	25	20	5	2.5
Deep Type 40mm	40	30	20	5	5
Deep Type 50mm	50	40	20	5	5
Deep Type 60mm	60	50	20	5	5
Deep Type 70mm	70	60	20	5	5
Deep Type 80mm	80	70	20	5	5
Deep Type 90mm	90	80	20	5	5
Deep Type 100mm	100	90	20	5	5
Deep Type 120mm	120	110	20	5	5